# **Universal Radiographic System**

### **SPECIFICATIONS**

#### Motorized swivel arm:

**Dual-speed** motorized movement for:

1) SID adjustment; 2) arm elevating movement

3) Rotation of the arm

Source-to-image disance (SID)

100-180 cm (39.4" - 70.9") continuosly adjustable

(dual-speed motorized movement)

Vertical travel Rotation angle 1.215 mm (47.8") 150° (90° +/- 30°)

#### Image detector cabinet:

Detector cabinet with integrated Control Panel - image digital detector: maximun size: 43 cm x 43 cm (17" x 17")

- range of tilt: +/- 45º

Grid type

1) Reciprocating bucky (oscillating grid) (103 lines / Ratio 12:1 / f = 150 cm / 59")

- 2) Grid cabinet + non oscillating interchangeable grid a) 1st, grid - 180 lines / Ratio 12:1 / f = 100 cm / 39.4" b) 2nd grid - 180 lines / Ratio 12:1 / f = 180 cm / 70.9"

#### Mobile patient table:

Table-top composition Patient capacity Dimensions (L x W x H) Laminated structure 200 kg (441 lb)

2.000 mm x 650 mm x 700 mm (78.7" x 25.6" x 27.6")

#### Options:

Mobile table Collimator Unit fixation system Image Detector Cabinet Carbon fibre table top - Patient capacity: 200 kg (441)

Automatic collimator Floor to celling support Motorized rotation movement

#### Main column:

Floor to wall mounted support Dual-speed motorized swivel arm vertical movement for effortless and acurate setting

#### X-Ray generator:

SEDECAL Compact X-ray generator (up to 80 kW output power) available in three power versions: line powered / battery powered / capacitor assisted

#### Collimator:

- Manual type collimator
- Linear laser beam for patient positioning in combination with light field indicator
- Rails for accesories
- Measuring tape

#### Weights:

- SEDECAL X Plus LP net wight: 324 kg (714.3 lb)
- (mobile table and digital image detector are
- Mobile table net weight: a) standard table: 46 kg (101.4 lb) b) carbon fibre table: 30 kg (66.1 lb)

#### Power line requeriments:

(SEDECAL X Plus LP only) Line Voltage: 208; 230; 240 VAC (single phase) Power consumption: 1,5 kVA

Specificifications subject to change without prior notice

#### **SEDECAL**

C/ Pelaya, 9 - Polígono Industrial "Río de Janeiro" 28110 - Algete (Madrid-SPAIN) Tel.: (34) 91 - 628 0544 / 91 - 628 1592 Fax: (34) 91 - 628 0574 Email: marketing@sedecal.com

SEDECAL USA Inc

2920 N Arlington Heights Rd. Arlington Heights, Illinois 60004 Tel.: 1 (847) 394 - 6960 Fax: 1 (847) 394 - 6966











## **Universal Radiographic System**




The C-arm desing maintains constant alignment between the x-ray tube and image receptor, regardless of C-arm tilt positions or image receptor angle.

Its extraordinary flexibility makes the system ideal for all patients in standing, sitting or laying position, including those who are disabled or physically restricted.









The Sedecal X plus LP provides a variable source to image distance (SID) of 100 to 180 cm (40" -72") to accommodate a range of radiographics studies.

The x-ray tube assembly can be easly moved from 100 to 180 cm (40" -72") at any C-arm angle position, by using the operator hand control.








#### Simple operation automatic operation

The movements can be controlled from four different places:

- 1 From the Touch Screen Panel outside of the Room at the working station
- **2** From the Touch Buttons located in the X-ray Tube cover, around the Touch Screen Display Control
- 3 From the Infrared Remote Control
- 4 From the Control Panel at the Detector side

During all movements a Button will need to be continuously pressed (dead-man) until reaching the final position. At this moment a message tells the operator that the position has been successfully reached.



#### Automation & simple operation

The idea under this new design is that just by selecting one patient from the Work List which will provide the specific study (anatomical region and view) to be performed in the examination, the System will be ready to move to a specific pre-programmed position, open the collimator shutters to the specific pre-programmed collimation and have the specific x-ray parameters for this kind of study.

As a safety measure, once a final selection of the Anatomical View is selected from the work list, the system will require the operator to confirm the movement to the pre-programmed position and pre-programmed collimation by pressing the bottom continuously until this position is reached.








### Easy positioning better productivity

The above describe procedure will enormously simplify any system positioning, providing greater productivity and allowing that not very experienced operators will easily auto-position the system. Patient positioning will also be faster and simpler to the operator as an actual photo of patient positioning can be displayed in both the working station and the touch screen in the x-ray tube cover for any specific study.






### **Digital Detector Assembly**



Sedecal X plus LP can be equiped with any digital x-ray detector up to 43 x 43 cm. (17" x 17") size.



TWO different GRIDS available for studies at different Focal Distances.

### **Automatic and Safe Positioning**



The system's operator-friendly controls dual speed motorized movements, and the intelligent anti-collision system (two laser beams + pressure sensor + bumper detector) makes the patient positioning quick, easy and safe.







### **Universal Radiographic System**



